

AQP11 Antibody

Catalog No: #43227

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | AQP11 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antigen affinity purification. |
| Applications | IHC |
| Species Reactivity | Hu |
| Specificity | The antibody detects endogenous levels of total AQP11 protein. |
| Immunogen Type | peptide |
| Immunogen Description | Synthetic peptide of human AQP11 |
| Target Name | AQP11 |
| Other Names | AQPX1 |
| Accession No. | Swiss-Prot#: Q8NBQ7 Gene ID: 282679 |
| Uniprot | Q8NBQ7 |
| GeneID | 282679; |
| Concentration | 1.8mg/ml |
| Formulation | Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol. |
| Storage | Store at -20°C |

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #43227 at dilution 1/30,

Background

Aquaporins facilitate the transport of water and small neutral solutes across cell membranes.

Note: This product is for in vitro research use only