

Rabbit (NE) ELISA Kit

Catalog No: #EK17955



Package Size: #EK17955-1 48T #EK17955-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rabbit (NE) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Rabbit

Application Details

Detect Range: 1.2µg/L -42µg/L

Sample Type: Serum, plasma, Urine, tissue homogenates and other biological fluids.

Sample Volume: 50ul

Assay Time: 3-5h

Detection wavelength: 450nm

Note: This product is for in vitro research use only