

PPP1R15A Antibody

Catalog No: #43273

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PPP1R15A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total PPP1R15A protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human PPP1R15A
Target Name	PPP1R15A
Other Names	GADD34
Accession No.	Swiss-Prot#: O75807 Gene ID: 23645
Uniprot	O75807
GeneID	23645;
Concentration	2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human prostate cancer tissue using #43273 at dilution 1/30,

Background

This gene is a member of a group of genes whose transcript levels are increased following stressful growth arrest conditions and treatment with DNA-damaging agents. The induction of this gene by ionizing radiation occurs in certain cell lines regardless of p53 status, and its protein response is correlated with apoptosis following ionizing radiation.

Note: This product is for in vitro research use only