

APCDD1 Antibody

Catalog No: #43286

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	APCDD1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total APCDD1 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human APCDD1
Target Name	APCDD1
Other Names	HHS; HTS; B7323; HYPT1; DRAPC1; FP7019
Accession No.	Swiss-Prot#: Q8J025Gene ID: 147495
Uniprot	Q8J025
GeneID	147495;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human tonsil tissue using #43286 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #43286 at dilution 1/20.

Background

This locus encodes an inhibitor of the Wnt signaling pathway. Mutations at this locus have been associated with hereditary hypotrichosis simplex. Increased expression of this gene may also be associated with colorectal carcinogenesis.

Note: This product is for in vitro research use only