

# Plant SA ELISA Kit

Catalog No: #EK18420



Package Size: #EK18420-1 48T #EK18420-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Plant SA ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Plant

## Application Details

Detect Range:0.5mg/L -50mg/L

Sample Type: Serum, plasma, Urine,tissue homogenates and other biological fluids.

Sample Volume:50ul

Assay Time: 3-5h

Detection wavelength:450nm

Note: This product is for in vitro research use only