

Plant JA ELISA Kit

Catalog No: #EK18624



Package Size: #EK18624-1 48T #EK18624-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Plant JA ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Plant

Application Details

Detect Range:30pmol/L-1000pmol/L

Sample Type: Serum, plasma, Urine,tissue homogenates and other biological fluids.

Sample Volume:50ul

Assay Time: 3-5h

Detection wavelength:450nm

Note: This product is for in vitro research use only