

## DSTYK Antibody

Catalog No: #43325

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	DSTYK Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DSTYK protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human DSTYK
Target Name	DSTYK
Other Names	RIP5; RIPK5; CAKUT1; DustyPK; HDCMD38P
Accession No.	Swiss-Prot#: Q6XUX3Gene ID: 25778
Uniprot	Q6XUX3
GeneID	25778;
Calculated MW	105kd
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:200-1:1000

Immunohistochemistry: 1:10-1:50

## Images



Gel: 6%SDS-PAGE


Lysate: 40 µg

Lane 1-2: 823 cells, hepg2 cells

Primary antibody: 1/200 dilution

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 1 minute



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #43325 at dilution 1/25,

## Background

This gene encodes a dual serine/threonine and tyrosine protein kinase which is expressed in multiple tissues. It is thought to function as a regulator of cell death. Multiple transcript variants encoding different isoforms have been found for this gene.?

Note: This product is for in vitro research use only