

MED4 Antibody

Catalog No: #43340

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MED4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total MED4 protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein of human MED4
Target Name	MED4
Other Names	ARC36; VDRIP; DRIP36; TRAP36; HSPC126
Accession No.	Swiss-Prot#: Q9NPJ6 Gene ID: 29079
Uniprot	Q9NPJ6
GeneID	29079;
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:20-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human colorectal cancer tissue using #43340 at dilution 1/25.



Immunohistochemical analysis of paraffin-embedded Human breast cancer tissue using #43340 at dilution 1/25.

Background

This gene encodes a component of the Mediator complex. The Mediator complex interacts with DNA-binding gene-specific transcription factors to modulate transcription by RNA polymerase II. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.?

Note: This product is for in vitro research use only