

## ZNF23 Antibody

Catalog No: #43350

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF23 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF23 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human ZNF23
Target Name	ZNF23
Other Names	KOX16; ZNF359; ZNF612; Zfp612
Accession No.	Swiss-Prot#: P17027 Gene ID: 7571
Uniprot	P17027
GeneID	7571;
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:100-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human Thyroid cancer tissue using #43350 at dilution 1/200,



Immunohistochemical analysis of paraffin-embedded Human Gastric cancer tissue using #43350 at dilution 1/200,

## Background

May be involved in transcriptional regulation. May have a role in embryonic development.

Note: This product is for in vitro research use only