

UQCRH Antibody

Catalog No: #43353

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	UQCRH Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total UQCRH protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein of human UQCRH
Target Name	UQCRH
Other Names	QCR6; UQCR8
Accession No.	Swiss-Prot#: P07919Gene ID: 7388
Uniprot	P07919
GeneID	7388;
Calculated MW	11kd
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Immunohistochemistry: 1:25-1:100

Images



Gel: 12%SDS-PAGE
 Lysate: 40 µg
 Lane: Mouse heart tissue
 Primary antibody: 1/200 dilution
 Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
 Exposure time: 15 seconds

Background

Cytochrome b-c1 complex subunit 6, mitochondrial?is a?protein?that in humans is encoded by the?UQCRHgene. This is a component of the

ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is part of the mitochondrial respiratory chain. This protein may mediate formation of the complex between cytochromes c and c1.

Note: This product is for in vitro research use only