

BBS9 Antibody

Catalog No: #43360

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | BBS9 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antigen affinity purification. |
| Applications | IHC |
| Species Reactivity | Hu |
| Specificity | The antibody detects endogenous levels of total BBS9 protein. |
| Immunogen Description | Fusion protein of human BBS9 |
| Target Name | BBS9 |
| Other Names | B1; D1; C18; PTHB1 |
| Accession No. | Swiss-Prot#: Q3SYG4Gene ID: 27241 |
| Uniprot | Q3SYG4 |
| GeneID | 27241; |
| Concentration | 0.8mg/ml |
| Formulation | Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol. |
| Storage | Store at -20°C |


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colorectal tissue using #43360 at dilution 1/200,



Immunohistochemical analysis of paraffin-embedded Human Lung tissue using #43360 at dilution 1/200,

Background

This gene is downregulated by parathyroid hormone in osteoblastic cells, and therefore, is thought to be involved in parathyroid hormone action in bones. The exact function of this gene has not yet been determined. Alternatively spliced transcript variants encoding different isoforms have been identified.?

Note: This product is for in vitro research use only