

Goat Anti-Human IgG Secondary Antibody AF594 Conjugated

Catalog No: #L3027



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|---------------------|--|
| Product Name | Goat Anti-Human IgG Secondary Antibody AF594 Conjugated |
| Host Species | Goat |
| Clonality | Polyclonal |
| Species Reactivity | Hu |
| Specificity | IgG (H+L) |
| Conjugates | AF594 |
| Excitation Emission | Excitation: 591nm Emission: 614nm |
| Formulation | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol |
| Storage | Store at -20 °C/1 year |

Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:800

Note: This product is for in vitro research use only