

# Goat Anti-Human IgG Secondary Antibody APC Conjugated



Catalog No: #L30216

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat Anti-Human IgG Secondary Antibody APC Conjugated
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	IgG (H+L)
Conjugates	Allophycocyanin (APC)
Excitation Emission	Excitation: 650nm Emission: 660nm
Formulation	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide
Storage	Store at -20 °C/1 year

## Application Details

IHC(1/100 - 1/1000), IF (1/100 - 1/1000), FC (1/1000 - 1/4000), ELISA (Use at an assay dependent concentration)

Note: This product is for in vitro research use only