

# Goat Anti-Mouse IgG Secondary Antibody AF680 Conjugated

Catalog No: #L3039



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat Anti-Mouse IgG Secondary Antibody AF680 Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Ms
Specificity	IgG (H+L)
Conjugates	AF680
Excitation Emission	Excitation: 684nm Emission: 702nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

## Application Details

Antibody Format: Whole IgG

Western Blot: 1:20,000-100,000

Note: This product is for in vitro research use only