

# Goat Anti-Mouse IgG Secondary Antibody HRP Conjugated

Catalog No: #L3032



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat Anti-Mouse IgG Secondary Antibody HRP Conjugated
Brief Description	Secondary Antibody
Clonality	Polyclonal
Species Reactivity	Ms
Specificity	IgG (H+L)
Conjugates	Horseradish Peroxidase
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 50% glycerol
Storage	Store at -20 °C/1 year

## Application Details

Antibody Format: Whole IgG

IHC/ICC: 1:500 - 1:5,000 ELISA: 1:5,000 - 1:50,000 Western Blot: 1:10,000 - 1:50,000

Note: This product is for in vitro research use only