

# Goat Anti-Rabbit IgG Secondary Antibody Cy3 Conjugated

Catalog No: #L30111



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody Cy3 Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Rb
Specificity	IgG (H+L)
Conjugates	Cyanine Cy3
Excitation Emission	Excitation: 550nm Emission: 570nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

## Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:800

Note: This product is for in vitro research use only