

Goat Anti-Rabbit IgG Secondary Antibody FITC Conjugated

Catalog No: #L30113



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|---------------------|--|
| Product Name | Goat Anti-Rabbit IgG Secondary Antibody FITC Conjugated |
| Host Species | Goat |
| Clonality | Polyclonal |
| Species Reactivity | Rb |
| Specificity | IgG (H+L) |
| Conjugates | Fluorescein (FITC) |
| Excitation Emission | Excitation: 492nm Emission: 520nm |
| Formulation | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol |
| Storage | Store at -20 °C/1 year |

Application Details

Antibody Format: Whole IgG

most applications: 1:100 - 1:200

Note: This product is for in vitro research use only