

Goat Anti-Rabbit IgG Secondary Antibody Rhodamine (TRITC) Conjugated

Catalog No: #L30114



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Goat Anti-Rabbit IgG Secondary Antibody Rhodamine (TRITC) Conjugated
Host Species	Goat
Clonality	Polyclonal
Species Reactivity	Rb
Specificity	IgG (H+L)
Conjugates	Rhodamine (TRITC)
Excitation Emission	Excitation: 550nm Emission: 570nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 10mg/ml Bovine Serum Albumin, 0.05% Sodium Azide, 50% glycerol
Storage	Store at -20 °C/1 year

Application Details

Antibody Format:Whole IgG
most applications: 1:100 - 1:200

Note: This product is for in vitro research use only