

FAM96A Antibody

Catalog No: #43380

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FAM96A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FAM96A protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein of human FAM96A
Target Name	FAM96A
Other Names	FA96A; FAM96A; family with sequence similarity 96, member A; FLJ22875; Protein FAM96A
Accession No.	Swiss-Prot#: Q9H5X1 Gene ID: 84191
Uniprot	Q9H5X1
GeneID	84191;
Calculated MW	18kd
Concentration	1.1mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1:1000

Images



Gel: 12%SDS-PAGE
Lysate: 40 µg
Lane: Human fetal intestines tissue
Primary antibody: 1/200 dilution
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution
Exposure time: 30 seconds

Background

May play a role in chromosome segregation through establishment of sister chromatid cohesion.

Note: This product is for in vitro research use only