

ABHD14B Antibody

Catalog No: #43384

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ABHD14B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ABHD14B protein.
Immunogen Description	Full length fusion protein of human ABHD14B
Target Name	ABHD14B
Other Names	CIB; HEL-S-299
Accession No.	Swiss-Prot#: Q96IU4 Gene ID: 84836
Uniprot	Q96IU4
GeneID	84836;
Concentration	1.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colorectal cancer tissue using #43384 at dilution 1/200.



Immunohistochemical analysis of paraffin-embedded Human Tonsil tissue using #43384 at dilution 1/200.

Background

ABHD14B is a 210 amino acid protein that localizes to both the cytoplasm and the nucleus where it exists as two alternatively spliced isoforms. Has hydrolase activity towards p-nitrophenyl butyrate (in vitro). May activate transcription.

Note: This product is for in vitro research use only