

FKBP7 Antibody

Catalog No: #43387

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FKBP7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FKBP7 protein.
Immunogen Description	Fusion protein of human FKBP7
Target Name	FKBP7
Other Names	FKBP23; PPIase
Accession No.	Swiss-Prot#: Q9Y680Gene ID: 51661
Uniprot	Q9Y680
GeneID	51661;
Concentration	1.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #43387 at dilution 1/40,

Background

The protein encoded by this gene belongs to the FKBP-type peptidyl-prolyl cis/trans isomerase (PPIase) family. Members of this family exhibit PPIase activity and function as molecular chaperones. A similar protein in mouse is located in the endoplasmic reticulum and binds calcium.?

Note: This product is for in vitro research use only