

CHGB Antibody

Catalog No: #43395

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CHGB Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total CHGB protein.
Immunogen Description	Synthetic peptide of human CHGB
Target Name	CHGB
Other Names	SCG1
Accession No.	Swiss-Prot#: P05060Gene ID: 1114
Uniprot	P05060
GeneID	1114;
Concentration	2.6mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:20-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #43395 at dilution 1/40,

Background

This gene encodes a tyrosine-sulfated secretory protein abundant in peptidergic endocrine cells and neurons. This protein may serve as a precursor for regulatory peptides.?

Note: This product is for in vitro research use only