

FOXN2 Antibody

Catalog No: #43397

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FOXN2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total FOXN2 protein.
Immunogen Description	Synthetic peptide of human FOXN2
Target Name	FOXN2
Other Names	HTLF
Accession No.	Swiss-Prot#: P32314Gene ID: 3344
Uniprot	P32314
GeneID	3344;
Concentration	0.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:100-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Brain tissue using #43397 at dilution 1/200.



Immunohistochemical analysis of paraffin-embedded Human Liver tissue using #43397 at dilution 1/200.

Background

This gene encodes a forkhead domain binding protein and may function in the transcriptional regulation of the human T-cell leukemia virus long terminal repeat.

Note: This product is for in vitro research use only