

## PALB2 Antibody

Catalog No: #43402

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PALB2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total PALB2 protein.
Immunogen Description	Synthetic peptide of human PALB2
Target Name	PALB2
Other Names	FANCN; PNCA3
Accession No.	Swiss-Prot#: Q86YC2Gene ID: 79728
Uniprot	Q86YC2
GeneID	79728;
Concentration	1.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1:20-1:100

## Images



Immunohistochemical analysis of paraffin-embedded Human lung cancer tissue using #43402 at dilution 1/30,

## Background

This gene encodes a protein that may function in tumor suppression. This protein binds to and colocalizes with the breast cancer 2 early onset protein (BRCA2) in nuclear foci and likely permits the stable intranuclear localization and accumulation of BRCA2.

---

Note: This product is for in vitro research use only