

ZNF365 Antibody

Catalog No: #43441

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	ZNF365 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total ZNF365 protein.
Immunogen Description	Synthetic peptide of human ZNF365
Target Name	ZNF365
Other Names	UAN; Su48; ZNF365D
Accession No.	Swiss-Prot#: Q70YC4Gene ID: 22891
Uniprot	Q70YC5
GeneID	22891;
Calculated MW	24kd
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:2000

Immunohistochemistry: 1:20-1:100

Images



Gel: 12%SDS-PAGE


Lysate: 40 µg

Lane: Human liver cancer tissue

Primary antibody: 1/200 dilution

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 seconds



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #43441 at dilution 1/20,

Background

This gene encodes several isoforms which have different expression patterns and functions. Mutation in this gene is associated with uric acid nephrolithiasis (UAN). Alternatively spliced variants, encoding distinct proteins, have been identified.

Note: This product is for in vitro research use only