

DAP Antibody

Catalog No: #43505

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DAP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DAP protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein
Target Name	DAP
Accession No.	Swiss-Prot#: P51397NCBI Gene ID: 1611
Uniprot	P51397
GeneID	1611;
Concentration	0.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1: 20-100

Images



The image is immunohistochemistry of paraffin-embedded Human liver cancer tissue using DAP Antibody at dilution 1/35. (Original magnification: x200)

Background

This gene encodes a basic, proline-rich, 15-kD protein. The protein acts as a positive mediator of programmed cell death that is induced by interferon-gamma. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

Note: This product is for in vitro research use only