

UBE2T Antibody

Catalog No: #43604

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UBE2T Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total UBE2T protein.
Immunogen Type	protein
Immunogen Description	Full length fusion protein
Target Name	UBE2T
Other Names	FANCT; PIG50; HSPC150
Accession No.	Swiss-Prot#: Q9NPD8NCBI Gene ID: 29089
Uniprot	Q9NPD8
GeneID	29089;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1: 20-100

Images



The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using UBE2T Antibody at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)

Background

The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. Defects in this gene have been associated with Fanconi anemia of complementation group T. Two transcript variants encoding different isoforms have been found for this gene.?

Note: This product is for in vitro research use only