

UBR2 Antibody

Catalog No: #43611

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	UBR2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total UBR2 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human UBR2
Target Name	UBR2
Other Names	C6orf133; bA49A4.1; dJ242G1.1; dJ392M17.3
Accession No.	Swiss-Prot#: Q8IWW8NCBI Gene ID: 23304
Uniprot	Q8IWW8
GeneID	23304;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1: 20-100

Images



The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using UBR2 Antibody at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using UBR2 Antibody at dilution 1/30, on the right is treated with fusion protein. (Original magnification: x200)

Background

This gene encodes an E3 ubiquitin ligase of the N-end rule proteolytic pathway that targets proteins with destabilizing N-terminal residues for polyubiquitylation and proteasome-mediated degradation. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only