

SERPING1 Conjugated Antibody

Catalog No: #C37916



Package Size: #C37916-AF350 100ul #C37916-AF405 100ul #C37916-AF488 100ul
 #C37916-AF555 100ul #C37916-AF594 100ul #C37916-AF647 100ul
 #C37916-AF680 100ul #C37916-AF750 100ul #C37916-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	SERPING1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SERPING1 protein.
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human serpin peptidase inhibitor, clade G (C1 inhibitor), member 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	C1IN; C1NH; HAE1; HAE2; C1INH
Accession No.	Swiss-Prot#:P05155NCBI Gene ID:710NCBI mRNA#:NCBI Protein#:NP_001073920/Q86U17
Uniprot	P05155
GeneID	710;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	55
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Background

This gene encodes a highly glycosylated plasma protein involved in the regulation of the complement cascade. Its protein inhibits activated C1r and C1s of the first complement component and thus regulates complement activation. Deficiency of this protein is associated with hereditary angioneurotic oedema (HANE). Alternative splicing results in multiple transcript variants encoding the same isoform.

Note: This product is for in vitro research use only