

ITGA10 Conjugated Antibody

Catalog No: #C44824



Package Size: #C44824-AF350 100ul #C44824-AF405 100ul #C44824-AF488 100ul
 #C44824-AF555 100ul #C44824-AF594 100ul #C44824-AF647 100ul
 #C44824-AF680 100ul #C44824-AF750 100ul #C44824-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	ITGA10 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms
Specificity	ITGA10 antibody detects endogenous levels of total ITGA10
Immunogen Description	A synthesized peptide derived from human ITGA10
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Integrin alpha-10 antibody Integrin, alpha 10 antibody ITA10_HUMAN antibody ITGA10 antibody OTTHUMP00000015588 antibody PRO827 antibody
Accession No.	Swiss-Prot#:O75578 NCBI Gene ID:8515 NCBI mRNA#: NCBI Protein#:
Uniprot	O75578
GeneID	8515;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	128
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only