

ZNF787 Antibody

Catalog No: #43790

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF787 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total ZNF787 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human ZNF787
Target Name	ZNF787
Other Names	TIP20
Accession No.	Swiss-Prot#: Q6DD87NCBI Gene ID: 126208
Uniprot	Q6DD87
GeneID	126208;
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1: 30-150

Images



The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ZNF787 Antibody at dilution 1/45, on the right is treated with synthetic peptide. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ZNF787 Antibody at dilution 1/45, on the right is treated with synthetic peptide. (Original magnification: x200)

Background

May be involved in transcriptional regulation.

Note: This product is for in vitro research use only