

GPR48 Conjugated Antibody

Catalog No: #C44931



Package Size: #C44931-AF350 100ul #C44931-AF405 100ul #C44931-AF488 100ul
 #C44931-AF555 100ul #C44931-AF594 100ul #C44931-AF647 100ul
 #C44931-AF680 100ul #C44931-AF750 100ul #C44931-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	GPR48 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms Rt
Specificity	GPR48antibody detects endogenous levels of total?GPR48
Immunogen Description	A synthesized peptide derived from human GPR48
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BNMD17 antibody G protein coupled receptor 48 antibody G-protein coupled receptor 48 antibody GPR48 antibody Leucine-rich repeat-containing G protein-coupled receptor 4 antibody Leucine-rich repeat-containing G-protein coupled receptor 4 antibody LGR4 antibody LGR4_HUMAN antibody
Accession No.	Swiss-Prot#:Q9BXB1NCBI Gene ID:55366NCBI mRNA#:NCBI Protein#:
Uniprot	Q9BXB1
GeneID	55366;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	104
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°Cin dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only