

## GPR26 Conjugated Antibody

Catalog No: #C44967



Package Size: #C44967-AF350 100ul #C44967-AF405 100ul #C44967-AF488 100ul  
 #C44967-AF555 100ul #C44967-AF594 100ul #C44967-AF647 100ul  
 #C44967-AF680 100ul #C44967-AF750 100ul #C44967-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GPR26 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms Rt
Specificity	GPR26antibody detects endogenous levels of total?GPR26
Immunogen Description	A synthesized peptide derived from human GPR26
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AXOR49 antibody DKFZp761H2121 antibody FLJ33774 antibody  G protein-coupled receptor 26 antibody G protein-coupled receptor 26 antibody G protein-coupled receptor GPR26 antibody G-protein coupled receptor 26 antibody Gpr26 antibody GPR26_HUMAN antibody MGC138216 antibody
Accession No.	Swiss-Prot#:Q8NDV2NCBI Gene ID:2849NCBI mRNA#:NCBI Protein#:
Uniprot	Q8NDV2
GeneID	2849;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	38
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°Cin dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

---

Note: This product is for in vitro research use only