

STK4 Conjugated Antibody

Catalog No: #C45098



Package Size: #C45098-AF350 100ul #C45098-AF405 100ul #C45098-AF488 100ul
 #C45098-AF555 100ul #C45098-AF594 100ul #C45098-AF647 100ul
 #C45098-AF680 100ul #C45098-AF750 100ul #C45098-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	STK4 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	STK4 Antibody detects endogenous levels of total STK4
Immunogen Description	A synthesized peptide derived from human STK4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	STK4, Mammalian sterile 20-like 1, MST-1, MST2, KRS2, Serine/threonine kinase 4, TIAC, STE20-like kinase MST1, YSK3, Kinase responsive to stress 2
Accession No.	Swiss-Prot#:Q13043NCBI Gene ID:6789NCBI mRNA#:NCBI Protein#:
Uniprot	Q13043
GeneID	6789;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	56
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only