

PCYT1A Conjugated Antibody

Catalog No: #C45181



Package Size: #C45181-AF350 100ul #C45181-AF405 100ul #C45181-AF488 100ul
 #C45181-AF555 100ul #C45181-AF594 100ul #C45181-AF647 100ul
 #C45181-AF680 100ul #C45181-AF750 100ul #C45181-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	PCYT1A Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	PCYT1A Antibody detects endogenous levels of total PCYT1A
Immunogen Description	A synthesized peptide derived from human PCYT1A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PCYT1A, CCT Alpha, CCTA, CT, CTA, CCT-alpha, CT A, PCYT1, CCT A, CTPCT
Accession No.	Swiss-Prot#:P49585NCBI Gene ID:5130NCBI mRNA#:NCBI Protein#:
Uniprot	P49585
GeneID	5130;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	42
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only