

BICD1 Antibody

Catalog No: #43846

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	BICD1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total BICD1 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human BICD1
Target Name	BICD1
Other Names	BICD
Accession No.	Swiss-Prot#: Q96G01NCBI Gene ID: 636
Uniprot	Q96G01
GeneID	636;
Concentration	0.4mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1: 20-100

Images



The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using BICD1 Antibody at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: x200)

Background

This gene is one of two human homologs of *Drosophila bicaudal-D*. It has been implicated in COPI-independent membrane transport from the Golgi apparatus to the endoplasmic reticulum. Two alternative splice variants have been described. Other alternative splice variants that encode different protein isoforms have been described but their full-length nature has not been determined.

Note: This product is for in vitro research use only