

KLHL36 Antibody

Catalog No: #43871

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KLHL36 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous levels of total KLHL36 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human KLHL36
Target Name	KLHL36
Other Names	C16orf44
Accession No.	Swiss-Prot#: Q8N4N3NCBI Gene ID: 79786
Uniprot	Q8N4N3
GeneID	79786;
Concentration	1.3mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1: 100-200

Images



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using KLHL36 Antibody at dilution 1/100, on the right is treated with synthetic peptide. (Original magnification: x200)



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using KLHL36 Antibody at dilution 1/100, on the right is treated with synthetic peptide. (Original magnification: x200)

Background

Probable substrate-specific adapter of an E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

Note: This product is for in vitro research use only