

GOLT1B Conjugated Antibody

Catalog No: #C45689



Package Size: #C45689-AF350 100ul #C45689-AF405 100ul #C45689-AF488 100ul
 #C45689-AF555 100ul #C45689-AF594 100ul #C45689-AF647 100ul
 #C45689-AF680 100ul #C45689-AF750 100ul #C45689-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	GOLT1B Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms
Specificity	GOLT1B Antibody detects endogenous levels of total GOLT1B
Immunogen Description	A synthesized peptide derived from human GOLT1B
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GOLT1B, Golgi transport 1B, GOT1A, HGOT1a, CGI-141, Golgi transport 1 homolog B, YMR292W
Accession No.	Swiss-Prot#:Q9Y3E0NCBI Gene ID:51026NCBI mRNA#:NCBI Protein#:
Uniprot	Q9Y3E0
GeneID	51026;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	15
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only