


## Description

Product Name	PPIG Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total PPIG protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human PPIG
Target Name	PPIG
Other Names	CYP; SRCyp; SCAF10; CARS-Cyp
Accession No.	Swiss-Prot#: Q13427NCBI Gene ID: 9360
Uniprot	Q13427
GeneID	9360;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: 1: 100-200

## Images



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using PPIG Antibody at dilution 1/100, on the right is treated with fusion protein. (Original magnification: x200)

## Background

PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides. May be implicated in the folding, transport, and assembly of proteins. May play an important role in the regulation of pre-mRNA splicing.

---

Note: This product is for in vitro research use only