

B3GALT5 Conjugated Antibody

Catalog No: #C45804



Package Size: #C45804-AF350 100ul #C45804-AF405 100ul #C45804-AF488 100ul
 #C45804-AF555 100ul #C45804-AF594 100ul #C45804-AF647 100ul
 #C45804-AF680 100ul #C45804-AF750 100ul #C45804-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	B3GALT5 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	B3GALT5 Antibody detects endogenous levels of total B3GALT5
Immunogen Description	A synthesized peptide derived from human B3GALT5
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	B3GALT5, B3T5, Beta-3-Gx-T5, Beta3GalT5, B3Gal-T5, B3galt5, B3GalT-V, Beta-1,3-GalTase 5, Beta3Gal-T5, GLCT5, B3GalTx
Accession No.	Swiss-Prot#:Q9Y2C3NCBI Gene ID:10317NCBI mRNA#:NCBI Protein#:
Uniprot	Q9Y2C3
GeneID	10317;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	36
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only