

DECR2 Antibody

Catalog No: #43885

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DECR2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total DECR2 protein.
Immunogen Type	protein
Immunogen Description	Fusion protein of human DECR2
Target Name	DECR2
Other Names	PDCR; SDR17C1
Accession No.	Swiss-Prot#: Q9NUI1NCBI Gene ID: 26063
Uniprot	Q9NUI1
GeneID	26063;
Calculated MW	31kd
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1000

Images



Gel: 8%SDS-PAGE
Lysate: 40 µg, Lane: HepG2 cell,
Primary antibody:DECR2 antibody at dilution 1/200 dilution,
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,
Exposure time: 10 seconds

Background

DECR2 (2,4-dienoyl-CoA reductase 2), also known as PDCR (peroxisomal 2,4-dienoyl-CoA reductase) or SDR17C1, is a 292 amino acid member of the short-chain dehydrogenases/reductases (SDR) protein family and the 2,4-dienoyl-CoA reductase protein subfamily. Localized to the peroxisome,

DECR2 is an auxiliary enzyme of beta-oxidation that catalyzes the NADP-dependent reduction of 2,4-dienoyl-CoA to yield trans-3-enoyl-CoA. DECR2 has also been shown to have catalytic activity towards 2,4,7,10,13,16,19-docosaheptaenoyl-CoA and short and medium chain 2,4-dienoyl-CoAs, suggesting that DECR2 is not a rate limiting step in the degradation of docosahexaenoic acid in the peroxisome. DECR2 is expressed as three isoforms produced by alternative splicing events.

Note: This product is for in vitro research use only