

CPD Conjugated Antibody

Catalog No: #C45850



Package Size: #C45850-AF350 100ul #C45850-AF405 100ul #C45850-AF488 100ul
 #C45850-AF555 100ul #C45850-AF594 100ul #C45850-AF647 100ul
 #C45850-AF680 100ul #C45850-AF750 100ul #C45850-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	CPD Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	CPD Antibody detects endogenous levels of total CPD
Immunogen Description	A synthesized peptide derived from human CPD
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CPD, Carboxypeptidase D, gp180, Metalloprotease D
Accession No.	Swiss-Prot#:O75976NCBI Gene ID:1362NCBI mRNA#:NCBI Protein#:
Uniprot	O75976
GeneID	1362;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	153
Formulation	Sodium Phosphate, 0.5% BSA, 40% glycerin, 0.02% sodium azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only