

PDCD6IP Antibody

Catalog No: #43906

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PDCD6IP Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total PDCD6IP protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human PDCD6IP
Target Name	PDCD6IP
Other Names	AIP1; ALIX; HP95; DRIP4
Accession No.	Swiss-Prot#: Q8WUM4NCBI Gene ID: 10015
Uniprot	Q8WUM4
GeneID	10015;
Calculated MW	96kd
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:200-1000

Images



Gel: 6%SDS-PAGE

Lysate: 40 µg, Lane: Hela cell lysate,

Primary antibody:PDCD6IP antibody at dilution 1/200 dilution,

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution,

Exposure time: 10 seconds

Background

This gene encodes a protein that functions within the ESCRT pathway in the abscission stage of cytokinesis, in intraluminal endosomal vesicle formation, and in enveloped virus budding. Studies using mouse cells have shown that overexpression of this protein can block apoptosis. In addition,

the product of this gene binds to the product of the PDCD6 gene, a protein required for apoptosis, in a calcium-dependent manner. This gene product also binds to endophilins, proteins that regulate membrane shape during endocytosis. Overexpression of this gene product and endophilins results in cytoplasmic vacuolization, which may be partly responsible for the protection against cell death. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. Related pseudogenes have been identified on chromosome 15.

Note: This product is for in vitro research use only