

P2RX1 Conjugated Antibody

Catalog No: #C46117



Package Size: #C46117-AF350 100ul #C46117-AF405 100ul #C46117-AF488 100ul
 #C46117-AF555 100ul #C46117-AF594 100ul #C46117-AF647 100ul
 #C46117-AF680 100ul #C46117-AF750 100ul #C46117-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	P2RX1 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms Rt
Specificity	P2RX1 Antibody detects endogenous levels of total P2RX1
Immunogen Description	A synthesized peptide derived from human P2RX1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	P2RX1, ATP receptor, p2X receptor 1, p2X1, p2X1 purinoceptor, p2X1 receptor, Purinergic receptor, p2RX1, p2X purinoceptor 1, p2X receptor, subunit 1
Accession No.	Swiss-Prot#:P51575NCBI Gene ID:5023NCBI mRNA#:NCBI Protein#:
Uniprot	P51575
GeneID	5023;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	45
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only