

# PIGQ Conjugated Antibody

Catalog No: #C46128



Package Size: #C46128-AF350 100ul #C46128-AF405 100ul #C46128-AF488 100ul  
 #C46128-AF555 100ul #C46128-AF594 100ul #C46128-AF647 100ul  
 #C46128-AF680 100ul #C46128-AF750 100ul #C46128-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PIGQ Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms Rt
Specificity	PIGQ Antibody detects endogenous levels of total PIGQ
Immunogen Description	A synthesized peptide derived from human PIGQ
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PIGQ, C407A10.1, PIG-Q, GPI1, HGPI1
Accession No.	Swiss-Prot#:Q9BRB3NCBI Gene ID:9091NCBI mRNA#:NCBI Protein#:
Uniprot	Q9BRB3
GeneID	9091;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	84
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

Note: This product is for in vitro research use only