

RCAN3 Conjugated Antibody

Catalog No: #C46652



Package Size: #C46652-AF350 100ul #C46652-AF405 100ul #C46652-AF488 100ul
 #C46652-AF555 100ul #C46652-AF594 100ul #C46652-AF647 100ul
 #C46652-AF680 100ul #C46652-AF750 100ul #C46652-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	RCAN3 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total RCAN3 protein.
Immunogen Description	Synthetic peptide corresponding to residues near the N terminal of human RCAN3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	RCN3; MCIP3; hRCN3; DSCR1L2
Accession No.	Swiss-Prot#:Q9JKK0NCBI Gene ID:53902NCBI mRNA#:NCBI Protein#:NP_038469
Uniprot	Q9JKK0
GeneID	53902;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	27
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

RCAN3 inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. It could play a role during central nervous system development.

Note: This product is for in vitro research use only