

Fascin 1 antibody

Catalog No: #22918

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	Fascin 1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 162 and 356 of Human FSCN1
Target Name	Fascin 1
Accession No.	Swiss-Prot:Q16658Gene ID:6624
Uniprot	Q16658
GeneID	6624;
Concentration	0.2mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 55kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
 A: HeLa
 B: Hep G2
 7.5% SDS PAGE
 Primary antibody diluted at 1: 1000

Background

Organizes filamentous actin into bundles with a minimum of 4.1:1 actin/fascin ratio. Probably involved in the assembly of actin filament bundles present in microspikes, membrane ruffles, and stress fibers.

Note: This product is for in vitro research use only