

## IFT22 Conjugated Antibody

Catalog No: #C47736



Package Size: #C47736-AF350 100ul #C47736-AF405 100ul #C47736-AF488 100ul  
 #C47736-AF555 100ul #C47736-AF594 100ul #C47736-AF647 100ul  
 #C47736-AF680 100ul #C47736-AF750 100ul #C47736-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IFT22 Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu, Ms, Rat
Specificity	The antibody detects endogenous levels of total IFT22 protein.
Immunogen Description	Fusion protein of human IFT22
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	RABL5
Accession No.	Swiss-Prot#:Q9H7X7NCBI Gene ID:64792NCBI mRNA#:NCBI Protein#:BC009823
Uniprot	Q9H7X7
GeneID	64792;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	21 kDa
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Background

---

Small GTPase-like component of the intraflagellar transport (IFT) complex B.

---

Note: This product is for in vitro research use only