

DYNAP Antibody

Catalog No: #44001

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | DYNAP Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antigen affinity purification |
| Applications | IHC |
| Species Reactivity | Hu |
| Specificity | The antibody detects endogenous levels of total DYNAP protein. |
| Immunogen Type | peptide |
| Immunogen Description | Synthetic peptide of human DYNAP |
| Target Name | DYNAP |
| Other Names | C18orf26 |
| Accession No. | Swiss-Prot#: Q8N1N2NCBI Gene ID: 284254 |
| Uniprot | Q8N1N2 |
| GeneID | 284254; |
| Concentration | 1.1mg/ml |
| Formulation | Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol. |
| Storage | Store at -20°C |

Application Details

Immunohistochemistry: 1: 30-150

Images



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using DYNAP Antibody at dilution 1/50, on the right is treated with synthetic peptide. (Original magnification: x200)

Background

Plays a role in the regulation of cell proliferation. Promotes activation of the AKT1 signaling pathway. Promotes phosphorylation of AKT1 at 'Ser-473'.

Note: This product is for in vitro research use only