

Melanocortin receptor 4

Catalog No: #AP79716



Package Size: #AP79716-1 50ug #AP79716-2 100ug #AP79716-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Melanocortin receptor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MC4-R
Accession No.	Uniprot ID:P70596
Uniprot	P70596
GeneID	25635;
SDS-PAGE MW	12.43kDa
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Expression Region:1-113aaStorage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mc4r

Note: This product is for in vitro research use only